GUTTA MASTER®

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Indication for Use

The GuttaMaster® Oven has been specially developed for heating endodontic obturators for use in root canal therapy. Do not use this oven for any other purpose.

Contraindications

There are none known.

Instructions for Use

1. Plug the oven into a 110/220 VAC outlet. Then turn the oven on with the power on/off switch located on the back. The yellow stand-by indicator will light up.
2. Put both obturator holders in the upper position by pressing the back of the $\hat{\uparrow}$ holder, and guiding the holder up with your finger.
3. After disinfecting the obturator and setting the rubber stop, place the obturator in the left obturator holder. **Regardless of what working length you have determined, the rubber stop must be under the holder.** Make sure the obturator is hanging straight, not at an angle. The obturator must clear all sides of the heating chamber.
4. Push the $\hat{\downarrow}$ holder down until you hear it click.
5. Push the button which corresponds to the size of the obturator you wish to heat. A yellow LED lights up. Then push the \odot LEFT button. The obturator is now being heated. The yellow LED is blinking. The heating time varies – depending on obturator size – from 20 to 45 seconds, and is regulated automatically.

The heating times for each obturator size are as follows:


Size	Signal Time	Tolerance
≤ 25	20 seconds	± 5 seconds
30-60	41 seconds	± 7 seconds
≥ 70	45 seconds	± 8 seconds

6. As soon as the obturator has reached the optimum temperature, a signal beeps and the yellow LED glows permanently. After the first signal “beep”, the obturator is ready for use. Push the $\hat{\uparrow}$ holder and guide it with your finger. Take the obturator carefully out of the holder by pulling it towards you, making sure not to scrape the obturator on any part of the holder.
7. You may leave the obturator in the GuttaMaster® Oven for up to 90 seconds after the first signal “beep”. It will keep the obturator at the appropriate temperature and ready for use. The oven will “beep” every 15 seconds to remind you that the obturator is still in the oven. After 90 seconds, the heating element will switch off automatically.

8. As an extra safety precaution, after you have taken the obturator out of the holder, switch off the heating element by pushing the  LEFT button again.
NOTE: Keep in mind that you won't be able to switch the heating element off during the heating cycle.
9. If you want to heat more obturators, you may alternate using the left and the right holders to continue your work efficiently. Wait until after the signal "beeps" for the first holder, then you may immediately start to heat the other one.
NOTE: It is not possible to heat two obturators at the same time. But you can leave the first obturator in the oven as long as the second obturator is heated up. Meanwhile, the obturator in the first holder will keep its optimum temperature for 90 seconds.
10. When you are finished using the oven, use the power on/off switch located at the back of the oven to turn it off. Then unplug the oven until you are ready to use it again.

Cleaning Instructions

Clean oven after each use. Any gutta-percha remaining on the holder or in the heating chamber may be effectively removed with Orange Solvent (generic dental solvent) on a cotton pellet.

Warning

1. Turn off and unplug oven when not in use.
2. Do not remove the cover of the oven.
3. This oven is designed for use in a dental office. Do not use in any other environment.

Precaution

1. To prevent fire or electric shock, do not expose this product to moisture.
2. To reduce the risk of electric shock, do not remove cover; there are no user-serviceable parts inside.
3. This oven is intended for use in heating gutta-percha coated obturators only. Do not insert anything other than an obturator into the oven.

Adverse Reactions

There are non known.

Technical Specifications

Voltage rating: 110/220 VAC, 50/60 Hertz, range of the supply voltage +/- 10%,
Fuse 2 Amps Slow.

Explanation of Marking Symbols

Alternating current

Power On


Power Off

Caution! Hot Surface

Do not dispose of as normal household waste;
special disposal according to the Directive 2002/96/EEC
of 23.01.2003 concerning waste of electrical/electronic equipment.

VW000264 Rev. 0/23.01.09

 **VDW GmbH**

P.O. Box 830954 • D-81709 München
Tel. +49 89 62734-0 • Fax +49 89 62734-304
www.vdw-dental.com • info@vdw-dental.com

Endo Easy Efficient®