DT ILLUSION™ SL

XRO® 


**SAFETY
LOCK**
TECHNOLOGY 

 **VDW**®

Endo **Easy** **Efficient**®

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INDICATIONS*

In case of insufficient residual tooth substance (more than one missing coronal wall), the post is needed to support the coronal restoration.

* Sorrentino R, Monticelli F, Goracci C, Zarone F, Tay FR, Garcia-Godoy F, Ferrari M
Effect of post-retained composite restorations on the fracture resistance of endodontically-treated teeth.
American Journal of Dentistry, 2007, 20: 269-274


Abb 1 / Fig. 1/ Schéma n°1/ rys 1/ obr 1/ рис 1 Abb. 2/ Fig. 2/ Schéma n°2/ rys.2/ obr.2/ рис. 2

DT ILLUSION™ SL

XRO® 

de	GEBRAUCHSANWEISUNG	1
en	INSTRUCTIONS FOR USE	10
fr	MODE D'EMPLOI	19
pl	INSTRUKCJA UŻYTKOWANIA	28
cs	NÁVOD K POUŽITÍ	37
ru	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ	46

Additional languages available upon request.

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Caution: For dental use only.

DESCRIPTION

DT ILLUSION™XRO®SL is a highly radiopaque (360 % - 730 % Al), translucent quartz fibre post in an epoxy resin-matrix. Its double taper permits a good fitting to the canal form as presented after modern treatment. This post design allows for inserting the post without unnecessary loss of tooth substance while offering optimal retention for building up the residual tooth.

Colour pigments with a thermal transition point are contained within the epoxy resin matrix (ILLUSION). They make the post change colour at ambient temperatures (22 °C) according to its colour coding and, for optimal aesthetics, the post becomes translucent at body temperature (37 °C) – the colour changes at 29 °C; cool air or water make the colour appear when necessary. It is possible to determine the size of the post simply on the basis of the colour. The colour provides improved visibility of the post during the treatment and makes it possible to clearly distinguish the post from the surrounding dentine. This is particularly helpful with retreatments*, as unnecessary loss of tooth substance or even a via falsa can be avoided when drilling the post. In addition, it makes it easier to see any possible remaining portions of the post.

To ensure the perfect aesthetics of the posts, the translucency transition point already starts at 29 °C, thus below body temperature (37 °C). The low transition point guarantees perfect aesthetics at all times.

*The use of DT Post Removal Kits is recommended for retreatments.

The surface of the DT ILLUSION™ XRO®SL quartz fibre post is pre-conditioned with silicate (SiOx) and silane (SiH4). A protective, laquerlike polymer layer ensures that the silicate and the silane are not contaminated or deactivated. This polymer layer is chemically integrated when inserting the post in the composite cement. The Safety Lock® coating does away with pretreating the post with silane or primer and ensures a long-term bonding of the post with the composite cement.

The initial length of the post is 20 mm. The apical length of the coated surface is 15 mm. For post dimensions refer to fig.2.

INDICATIONS FOR USE

In case of insufficient residual tooth substance (more than one missing coronal wall), the post is needed to support the coronal restoration. A ferrule effect (1.5 mm - 2 mm of coronal wall) should be possible. If not, a crown extending measure should be performed (surgically or orthodontically) in order to create a sound coronal dental structure.

CONTRAINDICATIONS

Tooth with insufficient residual circular coronal substance of 1.5 - 2 mm.

WARNING

None known.

PRECAUTIONS

- The post must be cleaned with alcohol before insertion into the canal.
- Do not touch the posts with your fingers after cleaning.
- **Do not** coat the post with silane.
- Shorten the post outside of the mouth.
- VDW recommends the use of a rubber dam whenever possible.

ADVERSE REACTIONS

According to our current standard of knowledge concerning that technique, none adverse reactions are known.

STEP-BY-STEP INSTRUCTIONS

1. Preparation, cleaning and obturation of the root canal
Prepare, clean, dry and fill the root canal according to the scientifically generally recognised techniques of endodontics.
Avoid the use of root canal filling materials containing eugenol. They can impair a hardening of the luting composites.
VDW recommends the use of a rubber dam whenever possible.
2. Select the correct DT ILLUSION™XRO®SL post according to the anatomical situation, the radiograph and the indications on fig. 1.

3. Determine the post length considering that approximately 1/2 to 1/3 of the total post length must be inserted into the root canal. Remove the root filling with the DT Universal Drill, ref. 0935 000 (rotation speed 1000-2000 rev/min.). At least 4 mm of root canal filling should remain apically. However, in the presence of severely curved canals, this calculated length cannot always be achieved. The unnecessary loss of tooth substance is to be prevented.
4. Shape the canal with the DT Finishing Drill (ref. 0935 001 or 0935 002 or 0935 003) according to the selected DT ILLUSION™ XRO®SL at a rotation speed between 1000-2000 rev/min. If post size 0 is used, no additional drilling should be necessary. The DT Universal Drill is the precision drill for size 0.
It should be ensured that all of the gutta-percha is removed from the canal walls to facilitate an adhesive bonding – this applies in particular to oval canals.
5. Rinse the root canal (aqua dest.) to remove dentine chips and dry with paper points.
6. Clean the post with alcohol before fitting.
7. Check if the post properly fits in the canal.
8. Shorten the coronal post to its final length outside of the mouth with a diamond disc. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
9. Clean the post with alcohol. Remove excess alcohol or allow it vaporise.

Adhesive bonding

10.DT ILLUSION™XRO®SL quartz fibre posts can be bonded with all of the customary adhesive procedures and all composite cements on a BisGMA- or UDMA basis. A total etch – total bond procedure is recommended for the highest bonding values.

The following products were successfully tested with the Safety Lock® coating:

- Adhesive systems: self- or dual curing:
Prime&Bond NT, XP Bond + SCA (Dentsply); Clearfil SE BOND (Kuraray); Syntac, ExciTE DSC, Adhese (Ivoclar Vivadent); Scotchbond (3M Espe); Gluma Comfort Bond (Heraeus Kulzer); ALL-BOND 2, ONE-STEP, ONE-STEP PLUS (Bisco)
- Composite cement: self- or dual curing:
Core X-flow, SmartCem 2, Calibra, X-Flow (Dentsply); Panavia F 2.0 (Kuraray); Variolink II, Multilink Automix (Ivoclar Vivadent); LuxaCore (DMG); Rely X Unicem, RelyX ARC (3M Espe); ParaCem (Coltene Whaledent); 2Bond2 (Heraeus Kulzer); Duo Link (Bisco); G-Cem (GC)
- Core build-up: composite material

Adhesive bonding using the example of the total etch – total bond procedure:

11. Apply etching gel* in the canal and on the exposed dentine in accordance with the manufacturer's instructions. Next rinse and dry gently with paper points, but leave canal wall surface moist. Avoid contact with gums – or, change for a self-etching primer.
12. Apply adhesive* (dual- or self-curing) with a very thin brush (e.g. VDW Composibrush) in the root canal and on the coronal dentine. Remove excess with paper points and gently air dry. The surfaces should have a uniform glossy appearance. If not, repeat application and drying as described above.
DT ILLUSION™ XRO®SL has a Safety Lock® coating, which is made of a silane and silicate layer, protected against contamination and deactivation by a polymer layer. The Safety Lock® coating does away with pretreating the post with silane or primer and ensures a long-term bonding of the post with the composite cement.
13. Mix the dual curing composite cement* (do not use any exclusively light-curing composite cements) and apply it in the canal using a tube, then apply composite cement on the post. Seat the post immediately. Remove excess with an appropriate instrument or a foam pellet. Afterwards light-cure and in the meantime stabilise the post by applying gentle pressure with the tip of the curing light fibre optic.
14. For a subsequent composite build-up of the exposed post, apply the bonding agent* with the post (only if less than 5 mm of the quartz fibre post has been shortened) and create the build-up of the post with the composite material. A one-stage procedure is likewise

* Please refer to the manufacturer's instructions for use.

possible (post insertion and post build-up in one go). It is recommended that the coronal end of the post is completely covered with the composite in order to avoid the possible penetration of liquid along the quartz fibre and thereby a hydrolytic degradation of the bonding agent. In cases where a prosthetic crown is not indicated, a minimum of 1 mm of composite restorative should cover the post. The composite should be appropriately-shaded, polishable and aesthetic.

Retreatment:

If a retreatment of the post should be necessary, it is recommended to use the DT Post Removal Kit (ref. 0938 001). The DT ILLUSION™XRO®SL quartz fibre posts can be drilled using the help of the DT Pilot Drill and the DT Carbide Drill.*

DT ILLUSION™XRO®SL offers a unique thermal colour adaptation, which enables better visibility of the post during retreatments and makes it distinguishable from the surrounding dentine. Thermally reactive colour pigments change the colour of the post at ambient temperature (22 °C) according to its colour coding. For optimal aesthetics it becomes translucent at body temperature (37 °C) – the colour changes at 29 °C; cool air or water make the colour appear if necessary. This is particularly helpful with retreatments, as unnecessary loss of tooth substance or even a *via falsa*, particularly when drilling the post, can thus be avoided. In addition, it makes it easier to see any possible remaining portions of the post.

* Observe the instructions for use of the DT Post Removal Kit.

STERILISING, DISINFECTING AND MAINTENANCE

Disinfect DT ILLUSION™XRO®SL with medical isopropyl alcohol.

Do not use a thermo disinfectant.

DT Universal Drill and DT Finishing Drill

The instruments are supplied unsterile. Please disinfect and sterilize instruments before each use. Disinfect instruments with disinfecting or cleaning agents appropriate for rotary instruments and follow the manufacturer's instructions for use. Make sure instruments do not touch in the ultrasonic cleaner. Do not disinfect instruments in thermo disinfectant. Sterilize instruments in autoclave (134 °C – 18 minutes – 1 cycle) or hot air sterilizer (180 °C).

Please refer to the device manufacturer's instructions for use. For better results, use autoclave bags or sterilizing bags for each drill. Do not store instruments in wet or humid conditions. Control instruments visually before each use: change instruments when they are damaged, corroded or dull. It is advised not to use the drill more than 15 times.

STORAGE

Keep out of direct sunlight. Store at room temperature.

SAFETY AND LIABILITY

VDW declines any liability or compensation for possible damage due to:

- use of instruments which do not belong to the system and which may impair the function.
- non-observation of the instructions for use.

The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the instructions for use.

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VDW GmbH
VDW GmbH • Bayerwaldstr. 15 • 81737 Munich • Germany