
DT Light[®] Post

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INDICATIONS*

In case of insufficient residual tooth substance (more than one missing coronal wall), the post is needed to support the coronal restoration.

- * Sorrentino R, Monticelli F, Goracci C, Zarone F, Tay FR, Garcia-Godoy F, Ferrari M
Effect of post-retained composite restorations on the fracture resistance of endodontically-treated teeth.
American Journal of Dentistry, 2007, 20: 269-274



Abb. 1 / Fig. 1/ Schéma n°1/ Rys.1 / Obr.1

Abb. 2 / Fig. 2 / Schéma n°2/ Rys.2 /Obr.2

DT Light[®]**Post**

Translucent quartz fibre technology

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DESCRIPTION

DT Light®Post is a radiopaque translucent quartz fibre post with a double tapered shape. Initial length of the post is 20 mm. For post dimensions refer to fig. 2.

INDICATIONS FOR USE

In case of insufficient residual tooth substance (more than one missing coronal wall), the post is needed to support the coronal restoration. A ferrule effect (on at least 2 mm of coronal wall) should be possible. If not, the tooth should be treated periodontically or orthodontically in order to create a sound coronal dental structure.

CONTRAINDICATIONS

None known to date.

PRECAUTIONS

- The post must be cleaned with alcohol before insertion into the canal.
- Do not touch the posts with your fingers after cleaning.
- Shorten the post outside of the mouth.
- Remove the coloured ring from the post before trying it in the canal.

ADVERSE REACTIONS

According to our current standard of knowledge none known.

STEP-BY-STEP INSTRUCTIONS

1. Select the correct DT Light®Post according to the anatomical situation, the radiograph and the indications on fig. 1. **Remove the coloured ring.**
2. Determine the post length considering that approximately 1/2 to 1/3 of the total post length must be inserted into the root canal. Remove the root filling with the DT Universal Drill, ref. 0935 000 (rotation speed 1000-2000 rpm). At least 4 mm of root canal filling should remain apically. However, in the presence of severely curved canals this calculated length is not always correct.
3. Shape the canal with the DT Finishing Drill (ref. 0935 001 or 0935 002 or 0935 003) according to the selected DT Light®Post at a rotation speed between 1000-2000 rpm. If post size 0 is used, no additional drilling should be necessary. The DT Universal Drill is the precision drill for size 0.
4. Check if the post is properly seated in the canal.

5. Shorten the post to its final length outside of the mouth with a diamond disc. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
6. Clean the post with alcohol. Remove excess.
7. Apply etchant* in the canal and on the exposed dentin and leave it for 15 seconds. Rinse for 10 seconds. Dry gently with paper points, but leave canal wall surface moist. Avoid contact with gums – or, change for a self-etching primer.
8. Mix the primer* and its self-cure activator* (the primer can be dual- or self-curing). With a very thin brush (e.g. VDW Composibrush) apply 2 coats of this preparation in the root canal and leave them for 20 seconds. **Remove excess with paper points** and gently air-dry all surfaces. The surfaces should have a uniform glossy appearance. If not, repeat application.
9. Apply a single coat of primer to the post. Gently air dry 5 seconds and light-cure the post for 10 to 20 seconds outside the mouth.

* Please refer to the manufacturer's instructions for use.

10. Mix the cement* (dual- or self-curing) and apply it in the canal and – if wished - on the post using a Lentulo spiral filler (e.g. VDW) or a syringe tip. Seat the post immediately. Remove excess with an appropriate instrument. Light-cure for 40 to 60 seconds, applying gentle pressure with the tip of the curing light fibre optic.
11. Apply 2 coats of bonding agent on the exposed post, the cement and the coronal areas and leave them for 20 seconds. Remove excess with air jet. Light-cure for 20 seconds.
12. Directly model the build-up using a composite core build-up material*.

OTHER RECOMMENDATIONS

We recommend sterilization in autoclave (max 134°C). Use an autoclave bag for each post. DT Light® posts can be disinfected with medical alcohol. Do not use a thermal disinfectant.

DT Universal Drill and DT Finishing Drill: the instruments are supplied unsterile. Please disinfect and sterilize instruments before each use. Disinfect instruments with disinfecting or cleaning agents appropriate for rotary instruments and follow the manufacturer's

* Please refer to the manufacturer's instructions for use.

instructions for use. Make sure instruments do not touch in the ultrasonic cleaner. Do not disinfect instruments in thermo disinfectant. Sterilize instruments in autoclave (134°C) or hot air sterilizer (180°C). Please refer to the device manufacturer's instructions for use. For better results, use autoclave bags or sterilizing bags for each drill. Do not store instruments in wet or humid conditions. Control instruments visually before each use: change instruments when they are damaged, corroded or dull.

VDW recommends the use of a rubber dam whenever possible.

SAFETY AND LIABILITY

VDW declines any liability or compensation for possible damage due to:

- use of instruments which do not belong to the system and which may impair the function.
- Non-observation of the instructions for use

The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the instructions for use.

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VDW GmbH nepřebírá záruky za škody, které vzniknou použitím nesystémových nástrojů nebo jiným použitím.
VDW GmbH
VDW GmbH • Postfach 830954 • D-81709 München