
Motor related functions


- 1. POWER:** ON/OFF
- 2. CAL:** Calibration of contra-angle: at least once a day, after changing the contra-angle and after each lubrication.
- 3. ASR:** **Green:** automatic reverse motion as soon as torque limit has been reached.
Red: only reverse motion.
OFF (no LEDs on): automatic stop as soon as torque limit has been reached.
- 4. ANA:** **Green:** general reduction of torque values: recommended for difficult canal anatomy and for beginners.
- 5. MENU:**
 - Activation/deactivation of file systems.
(Upon delivery M_{two}®, FlexMaster®, and Dr's Choice are activated and thus visibly listed in the display, the other systems are invisible in the background. To activate/inactivate each file system, select it on the MENU and set it to ON or OFF. Only activated file systems are listed in the display.)
 - Activation/deactivation of the automatic apical stop when the apex has been reached (combined length determination with preparation).
 - Adjustment of sound level of motor and apex locator separately.
 - Selection of language: German or English.
- 6. ✓ CONFIRM:** For confirmation of all changes, e.g. torque and speed.
As alternative to the foot pedal:
Starts file rotation if the ✓ key is being pressed and held for 1.5 seconds.
Stops rotation if any other key has been pressed.
- 7./8. |◀/▶|:** Scrolls to the right/left: selection of system, file type, torque and speed.

9./10. +/- : Scrolls up and down:
Increase/reduction of torque and speed values.
Selection of next file system (Dr's Choice program as well) or next file type.

11. BATTERY: **Green:** battery fully charged
Red flashing: battery nearly empty – charging needed
Orange flashing: battery is recharging

Navigation in display

Upper row: shows the selected file system.

Lower row:

- a) sys** file system (e.g. M_{two}®, FlexMaster®, Dr's Choice)
- b) file** file type (10/04, 15/05...)
- c) gcm** torque in gcm (individually adaptable in all file systems)
- d) rpm** speed in rpm (individually adaptable in all file systems)


|◀/▶| keys: Jumps to the right/left between above shown display fields.

(▶ Arrow marks each active field which will be changed.)

+/- keys: Scrolls up/down through systems and file types to change torque and speed values.

Reset: POWER key and ✓ CONFIRM key simultaneously pressed and held: all individually changed values and activated file systems are reset to the values upon delivery.

Plugs (back of device)

Functions of apex locator

Important advice for successful length determination

- Use only the original contra-angle 6:1!
- For insulation purpose always wear protection gloves, use a rubber dam and/or the contra-angle sleeve!
- If the result is implausible, exclude probable causes (see "trouble shooting" in the user manual).
- In rare cases, the patient may feel pain: interrupt the treatment and place a medical dressing (see chapter 4.9 in the user manual).

• Connectors for cables (see front page)

Plug the lip clip cable properly in B and the file clamp cable in C. Check the connections!

• Switch-on of the apex locator (inclusive cable check)

Connect file clamp and lip clip

Stand-alone length determination with a hand file (without preparation)

or Connect contra-angle with lip clip

Combined length control with the contra-angle

• Stable connection

The display reads "apex locator ON" and the APEX LED switches to green.


• Length determination

- Create an initial gliding path with a hand instrument and determine the length in the stand-alone mode with VDW.GOLD® using a hand file and the lip clip.
- Dry the access cavity with a blower or a cotton pellet to avoid leakage current.
- Use the simultaneous length control with the contra-angle and the lip clip.
- If necessary, dry the canal, as too much liquid will cause erroneous results.

• Apical automatic stop


Upon delivery this auto-function is inactive. To activate this apical auto stop enter MENU and choose "Apex stop ON". When activated, the file tip stops automatically at the apex. By releasing and repressing the foot pedal, the rotation starts in reverse motion.

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