VDW.ULTRA®

Gebrauchsanweisung
Operating Manual
Manual de instrucciones
Manuale d'uso
Mode d'emploi

Endo Easy Efficient®


Abb. 1 Gesamtansicht – Overview – Vista general – Veduta d'insieme – Vue d'ensemble

	de	en	es	it	fr
1	Fußpedal	Foot pedal	Interruptor de pedal	Interruttore a pedale	Pédale
2	AN/AUS Anzeige	ON/OFF LED	Pulsador ON/OFF	Manopola di comando ON/OFF	Interrupteur ON/OFF
3	Einstellknopf	ON/OFF adjustment knob	Botón de regulación	Manopola di regolazione	Bouton de réglage
4	Handstückablage	Handpiece holder	Bandeja para la pieza de mano	Supporto per il manipolo	Support pour pièce à main
5	Handstück	Handpiece	Pieza de mano	Manipolo	Pièce à main
6	Handstückkabel	Handpiece cord	Cable de la pieza de mano	Cavo del manipolo	Cordon de pièce à main
7	Einstellknopf für Wasserzufuhr	Water spray adjustment knob	Botón de regulación para suministro de agua	Manopola di regolazione adduzione acqua	Bouton de réglage du débit d'eau


Abb. 2 Vorderansicht – Front panel – Vista frontal – Lato anteriore – Face avant


Abb. 3 Rückansicht – Rear panel – Vista trasera – Lato posteriore – Face arrière

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Contents

Introduction	27	APPENDIX	
1. Safety Information	28	11. Electromagnetic Compatibility	43
2. Standard Components	30	11.1. Electromagnetic Emissions	43
3. Description	30	11.2. Electromagnetic Immunity	44
4. Step-by-step Instructions	31	11.3. Electromagnetic Immunity/ Mobile Radiofrequency Equipment	45
4.1. Installation	31	11.4. Recommended Separation Distances	46
4.2. Before Operating the Unit	31	11.5. Cable Lengths	47
4.3. Operating the Unit	32		
4.3.1. Endodontic Applications	32		
4.3.2. Other Applications	34		
4.4. After Operating the Unit	34		
5. Sterilization and Replacement Parts	34		
5.1. Autoclavable Components	34		
5.1.1. Handpiece	34		
5.1.2. Tips	35		
5.2. Non-Autoclavable Components	35		
5.3. Replacement Parts	35		
6. Maintenance	36		
6.1. Handpiece Maintenance	36		
6.2. Tip Maintenance	36		
6.3. Control Unit Maintenance	36		
6.4. Water Filter Replacement and Maintenance	37		
7. Troubleshooting	38		
8. Technical/Performance Data	39		
9. Disposal and Recycling	40		
10. Warranty incl. Service Form	40		

Thank you for purchasing the VDW.ULTRA® ultrasonic device.

Indications for use

The VDW.ULTRA® is a piezo-electric ultrasonic generator that uses ultrasonic tips to carry out dental treatments, especially designed for root canal treatments.

In order to receive the most out of this high technology product, and prolong its life, we ask you to read these instructions carefully before using the device for the first time.

The User Manual is supplied for safety reasons and must be considered as an integral part of the system. We recommend always keeping it close at hand.

It is of the utmost importance that this manual be retained for future consultation. The manual must accompany the system in all cases of sale or other transfer in order that you as the new owner can be aware of its operational characteristics, the relative precautions and warnings and the terms of warranty.

Please do not hesitate to contact VDW GmbH for help with any doubt or problem that may arise during consultation of this User Manual.

See Directions

This symbol found on labels of the unit and accessories reminds you to consult this User Manual.

Notes

- *These instructions for use are available in several languages upon request.*
- *Changes to the content of these instructions are reserved without prior notice.*

Distributor:
VDW GmbH
Bayerwaldstr. 15
81737 Munich
Germany
www.vdw-dental.com
e-mail: info@vdw-dental.com

Manufacturer:
Satelec®, a company of ACTEON group
17 Ave. Gustave Eiffel, BP 30216
33708 Merignac Cedex
France
www.acteongroup.com

1. Safety Information

Note the meaning of the following symbols:

	Alternating current
	Protection class II
	Type BF device
	See Directions of Use; if the instructions are not being followed properly, operation may result in hazards for the product or the user/patient.
	Warning or Precaution
	Additional information, explanation on operation and performance.
	Do not dispose of as normal household waste; special disposal according to the Directive 2002/96/EEC of 27.01.2003 concerning waste of electrical/electronic equipment
	CE marking
	Autoclavable at 134° C
	Manufacturer's batch code
	No of items per package

Most operation and maintenance problems result from insufficient attention being paid to basic safety precautions and not being able to foresee the possibilities of accidents.

Problems and accidents are best avoided by foreseeing the possibility of danger and operating the unit in accordance with the manufacturer's recommendations.

Contraindication

Use of the VDW.ULTRA® is contraindicated in patients and/or the operator with an active implant such as pacemaker, etc.

General Precautions

Take the following precautions to reduce the risk of any accident:

- The VDW.ULTRA® is designed for dental use only by trained and qualified practitioners such as dentists and hygienists who act according to the national regulations.
- The device has been developed to be used with VDW endodontic tips.

For hygienic and prophylactic purposes Satelec® periodontic and scaling tips can also be used with the VDW.ULTRA®.

Using non VDW/Satelec® tips may damage the handpiece and even the tips themselves.

- Always place a rubber dam to avoid inhaling or swallowing any parts.
- The electrical supply to which the device is connected should comply with standards in force.
- The water supply system should satisfy the quality criteria required for the practice of dental treatment. It should be equipped with a supply cut-off system.
- Do not cover the device or obstruct the air vents.
- Do not use the device if it appears to be damaged or defective.
- Do not immerse in liquid or use outdoors.
- Do not place the device near a heat source.

- Make sure that the cords do not hinder the free passage of people.
- Before unplugging the power cord, turn-off the water supply and switch the device off.
- After use, turn-off the water supply.
- When unplugging the power cord, hold both the cord plug and the power socket in the wall.
- The device must not be used in the presence of explosive anaesthetic gases.
- The device should be stored in its original packaging, in an appropriate place, where it will not cause any danger to people.
- Do not carry out repairs or modifications to the device without prior authorization by VDW GmbH. If any fault occurs, contact your local supplier rather than have it repaired by an unauthorized person.
- The device must not be connected to or used in combination with any other apparatus or system. It must not be used as an integral component of any other apparatus or system. VDW GmbH will not be responsible for accidents, equipment damage, bodily injury or any other trouble which results from ignoring this prohibition.
- In case of any doubt, contact your local supplier or the After Sales Service of VDW GmbH in Munich.
- The device complies with the Electromagnetic compatibility standards (CEI 60601-1-2), however, make sure that any possible electromagnetic interference does not present any additional risk.

Adverse reactions

Unknown.

2. Standard Components

Before set-up, check that all parts are present in the box (see Fig. 1):

- Central unit with cords and foot pedal
- Handpiece
- Tip wrench
- User manual

Notes

- *None of the components are delivered disinfected or sterilized!*
- *Prior to installation, check the device for a damage. Report any damage sustained during shipment to your local supplier within 24 hours of receipt of the device.*

high to very high from 60 to 100:
abbreviation MAXI stands for maximum
power

Notes

- *See detailed directions for use of tips in section 4.3. "Operating the Unit".*
- *Specific information on VDW endo tips can be found in the VDW Endo Tip Guide.*
- *This device has been designed especially for endodontic root canal treatment, but can also be used for scaling and periodontic treatments. For details see section 4.3.2. "Other Applications."*

3. Description

Front panel (see Fig. 2)

The power adjustment knob is used for switching the unit ON/OFF (ON is indicated by the green LED) and setting the ultrasonic power levels from 10 to 100 according to the recommended power range for each tip.

The VDW.ULTRA® can be set to one of the following power ranges when used with the VDW endodontic tips.

low from 10 to max 30:
abbreviation IRRI stands for irrigation
activation

low to medium from 30 to max 50:
abbreviation REDO stands for
retreatment

medium from 40 to max 60:
abbreviation CAVI stands for
access cavity

Rear panel (see Fig. 3)

There are four connections and one knob on the rear panel.

- The power cord to connect the device to the electrical supply (11).
- The foot pedal cord to operate the unit (10).
- The water supply tube including a filter (9).
- The handpiece cord to connect the device to the handpiece (8).

On the right hand side there is the water flow adjustment knob (7).

4. Step-by-step Instructions

4.1. Installation

Notes

- *The VDW.ULTRA® unit should be connected to the electrical supply in accordance with current standards.*
- *The water supply system should satisfy the quality criteria required for dental treatment. In addition, it should be equipped with a water supply back flow prevention system.*

1. Carefully remove the unit and the accessories from their packing and place them properly on a flat surface in the desired location. Check that the power adjustment knob is set to position OFF.
2. Contact your local Service Technician in order to connect the water tube to the water supply.
3. Connect the power plug to an electrical socket.
4. Place the foot pedal properly to make operation easy.
5. Connect the handpiece carefully by plugging it into the connector cord, aligning the electrical contacts (see Fig. 4 below).

Warning

Neither rotate the handpiece nor twist the handpiece connector on its cord.


Fig. 4 Handpiece Connection

6. Place the handpiece with the head ahead on the handpiece holder (holder pivotable by 60°).

4.2. Before Operating the Unit

1. Sterilize the handpiece after disconnecting (see section 5 "Sterilization").
2. Make sure there is no trace of moisture in the handpiece connection, if this is the case then wipe it and dry it with the airblower. Then reconnect the handpiece carefully to the handpiece connection without twisting and rotating.
3. Sterilize the tips prior to first use (see section 5.1. "Sterilization"). Select the appropriate tip, and connect it by screwing it on with the supplied wrench. Tighten it moderately (see Fig. 5 below).


Fig. 5 Tip Tightening with Wrench

4. Check that all connections are properly made.
5. Prior to use, run a sufficient amount of water for several minutes to clear the line by pressing down the foot pedal with the power adjustment knob at the lowest level 10.

Warnings

- *Correct tip vibration requires a tight fit, but without over-tightening as there is no stop point. Over-tightening the tip with force may damage tip and/or handpiece.*
- *Do not rotate the handpiece during installation of the tip when the handpiece is already connected to the cord.*
- *Do not disconnect the handpiece when power is switched on and the foot pedal is pressed down.*

4.3. Operating the Unit

General Notes

- *The advanced technology of the device guarantees a stable power transfer to the tip to work at a constant amplitude.*
- *The indicated power ranges are only relevant for operating the VDW.ULTRA®. These recommended ranges may not be appropriate for other ultrasonic piezo-generators.*
- *The recommended settings are determined both by the clinician's experience with ultrasonics and the power needed to safely complete the specific procedure.*
- *The ultrasound tips are equipped with a water port to run water through the water supply tube in order to cool during the treatment.*

4.3.1. Endodontic Applications

Notes

- *Prior to the endodontic treatment place a rubber dam to avoid inhaling any foreign objects.*
- *Sufficient optical magnification is recommended.*

Precautions

- *Dry operation without cooling can lead to overheating of the handpiece and might cause a necrosis.*

To avoid this from happening, work intermittently, i.e. take breaks (2 mins operation/ 5 mins pause) when using the tip. For cooling external irrigation with a syringe is recommended.

- *For high-power applications, e.g. dislodging a metal post, it is mandatory to use water for cooling by adjusting the water spray knob.*

Other medium-power applications, e.g. during retreatment, can be carefully carried out in dry procedure for a short time.

Consider the following factors when selecting the endodontic tip:

1. Choose the tip size in accordance to coronal access, root anatomy, procedural application and general visibility.
2. Choose the appropriate minimum power setting for the tip that is being utilized.
3. Increase the power gradually. Observe whether this setting allows safe procedure. If not efficient enough, increase the power gradually up to the maximum power setting for that tip, as needed to maximize clinical performance.
4. Use irrigation whenever the clinical situation is appropriate.

Precautions regarding breakage

- *Do not switch quickly to the maximum level for each indication in order to avoid tip breakages. Too much vibration will damage the tip.*
- *Never activate an irrigation file by pressing the foot pedal as long as it is still in the air and has not yet got in touch with the root canal.*
- *Even if used according to the instructions – within the recommended power range – the tip can break in some cases.*

General procedure (but follow the detailed instructions in the VDW Endo Tip Guide):

Range IRR1

Activation of irrigation

1. Irrigate with an appropriate rinsing solution, e.g. a syringe with NaOCl or CHX. If the water spray was activated prior to rinsing, then close by adjusting the water spray knob.
2. Turn the power knob from low up to level 10.
3. Insert the tip into the root canal before activating.
4. Press the foot pedal for activating.
5. Check if power is sufficient. If not, gradually turn up.
6. We recommend activating of irrigation 3 times for around 20 seconds. Apply new rinsing solution each time.
7. Release the foot pedal before removing the tip from the canal.

 *Thin ultrasonic files can break – even if used according to the instructions.*

Range REDO

Removal of different materials during retreatment

Check whether cooling by adjusting the water spray knob is appropriate.

To avoid overheating, work intermittently by taking breaks.

1. Start at lowest level.
2. Press the foot pedal and check if power is sufficient.
3. If not, gradually turn up.

Range CAVI

Preparation of access cavity and pulp chamber floor

(e.g. troughing, removal of pulp stones)

Check whether cooling by adjusting the water spray knob is appropriate.

To avoid overheating, work intermittently by taking breaks.

1. Start at lowest level.
2. Press the foot pedal and check if power is sufficient.
3. If not, gradually turn up.

 *Not all CAVI tips can be used within the whole power range.*

Range MAXI

Only for dislodging metal posts with the MAXI tip (when a maximum tip vibration is indicated)

1. For cooling adjust the water spray knob.
2. To avoid overheating, work intermittently by taking breaks.
3. Start at lowest level.
4. Press the foot pedal and check if power is sufficient.
5. If not, gradually turn up.

Important Notes

- Always start at low level and slowly turn up!
- For details of the exact power range (e.g. Power 30-45) see the packaging label of each tip.

4.3.2. Other Applications

VDW.ULTRA® allows to carry out scaling and periodontic treatments with Satelec® tips which are compatible with the device.

In general, the scaling tips are used at high power, whereas the perio tips are used at low power.

To find the recommended settings on your VDW.ULTRA®, refer to the Equivalence Table below.

For example: the scaling tip no. 1 is recommended at a setting of 14-15 on the P5 Newtron ultrasound device from Satelec® which is equivalent to the setting of 70 to 75 on the VDW.ULTRA®.

VDW.ULTRA®	5	10	15	20	25	30	35	40	45	50
Satelec® P5 Newtron	1	2	3	4	5	6	7	8	9	10

VDW.ULTRA®	55	60	65	70	75	80	85	90	95	100
Satelec® P5 Newtron	11	12	13	14	15	16	17	18	19	20

Notes

- For further information contact Satelec®.
- Refer to the Tip Book of Satelec® which gives detailed instructions for proper settings and applications of scaling and perio tip.
- This warranty does not cover any damages caused by the application with Satelec® scaling and perio tips.

Precautions

- It is mandatory to always use the water spray during scaling and perio treatments.
- Do not switch quickly to the maximum level for each indication in order to avoid tip breakages. Too much vibration will damage the tip.

- After certain time of use the tips will change due to average wear. Refer to the tip card of Satelec®.

- Follow the instructions and precautions in the Tip Book from Satelec®.

4.4. After Operating the Unit

1. Turn-off the water supply
2. Switch-off the device by turning the knob to OFF.
3. Remove tip and disconnect handpiece.

5. Sterilization and Replacement Parts

Notes

- Follow the maintenance and/or sterilization instructions before using the accessories supplied for the first time.
- Avoid using cleaning agents and disinfectants containing inflammable substances. If unavoidable, ensure that the product has evaporated and that there is no combustible product on the device and accessories before operation.
- Consider concentration and contact time recommended by the manufacturer of the disinfectant.
- Use a disinfectant with tested effectiveness (e.g. DGHM certification or CE mark) or which is in compliance with any standard required by the national regulations.

Handpiece (except handpiece cord), wrench and tips are sterilizable.

5.1. Autoclavable Components

5.1.1. Handpiece

1. Flush the water circuits of handpiece and tip at low power level for around 20-30 sec. by depressing the foot pedal.
2. Disconnect both the tip and the cord.
3. Disassemble the nozzle from the handpiece.
4. Clean the handpiece with an alcohol wipe before sterilization.
5. Put it in a sterilization pouch for autoclaving.
6. Autoclave at 134°C for 18 minutes.
7. After sterilization thoroughly dry the electrical contacts of the handpiece.

5.1.2. Tips

Prior to sterilization carry out pre-disinfection, mechanical cleaning and drying.

1. Put on thick gloves when cleaning the tips with alcohol wipes or immerse the tips in an ultrasonic tank (enzymatic, quaternary ammonium base detergent/disinfectant solution).
2. Place under running water for at least 30 seconds.
3. Scrub the tips in a new disinfectant bath (enzymatic, quaternary ammonium base detergent/disinfectant solution) for at least 30 seconds, and in any case until all remaining traces of contamination have been removed.
4. Rinse again under running water for at least 30 seconds.
5. Finally, dry with clean, non-woven single-use cloth to remove any traces of liquid.
6. Put it in a sterilization pouch for autoclaving.
7. Autoclave at 134°C for 18 minutes.

Notes

- *Autoclave in a standard autoclave (type B; compliant with standard EN 13060) with a cycle of 134°C for 18 minutes at a pressure of 2 bar minimum.*
- *As not all autoclaves can reach 134°C and draw a pre-vacuum, please refer to your autoclave manufacturer for specific sterilization instructions.*
- *Autoclave and dry temperatures must not exceed 135°C.*
- *Thoroughly clean the components before autoclaving as described above.*
- *If chemical solutions or foreign debris are not removed, autoclaving could damage or deform the components.*
- *Autoclave handpiece and tips after each patient.*
- *A second handpiece can be purchased to ease the procedure.*
- *Do not sterilize in another way than in an autoclave.*
- *During sterilization, metal parts of different types must not touch each other. Any contact will lead to the formation of electrolytic coupling causing localized deterioration.*

5.2. Non-Autoclavable Components (Handpiece Cord and Control Unit)

Wipe handpiece cord and control unit with an appropriate dental disinfectant solution (e.g. Septol wipes).

5.3. Replacement Parts

- Replace the parts if necessary depending on degree of wear and duration of usage.
- Reorder replacement parts from your local supplier or directly at VDW GmbH in Munich.

6. Maintenance

6.1. Handpiece Maintenance

- As a routine: after treatment with water flow, flush the water circuits of the handpiece and tip with water at a low power level for around 20- 30 seconds by pressing down the foot pedal.
- In case of damage, excessive wear at least once a year: replace the watertight O-ring which is located in the rear part by using the fitting tool (for details see Fig. 6).

Note

Lubricate the water-tight O-ring with appropriate dental lubricant in order to prolong its efficacy and prevent internal leaks.


Fig. 6 Handpiece Maintenance

6.2. Tip Maintenance

Basically, shape and mass of a tip are the determining factors to obtain consistent performance of the ultrasonic generator.

- The thinner the diameter of a tip, the less the tip should be re-used.
- After a certain time of usage the shape of a tip will gradually be modified, i.e. the tip will become shorter.
- Operation with a worn tip will lead to longer preparation time and may increase temperature in the cavity or root canal. It may also lead to tip breakage.

1. Do not actively modify a tip by filing.
2. Regularly check a tip for visible wear with a loup or microscope.
3. Replace
 - the tip which is 3 to 4mm short after several operations, compared to the length of a new tip.
 - the tip which has apparently suffered an accidental impact (fall etc.).
 - the prebent tip (that has matched the canal curvature) already after 1 operation as prebending considerably modifies a tip and makes it a single use instrument.
 - the diamond coated tip as soon as it is blunt.

6.3. Control Unit Maintenance

1. Regularly inspect accessories and cords for insulation defects and replace if necessary.
2. Regularly check whether the air vents of the control unit are clean to avoid overheating.
3. It is necessary to remove any traces of moisture in the handpiece connection; if yes dry it by wiping and airblowing.

Warning

Do not use the device if it appears to be defective.

In case of any malfunction, contact the local supplier or directly VDW GmbH rather than have it repaired by an unauthorized person.

6.4. Water Filter Replacement and Maintenance

Regularly clean and change the water filter if necessary.

Follow the instructions below:

1. Turn-off the water supply and switch-off the device (OFF position), then unplug the power cord.
2. Unscrew the two parts of the filter using two open-ended 10 mm wrenches.


Fig. 7 Disassembly of Filter

3. Remove the filter cartridge to replace or clean it under running water.
4. Proceed in the same manner for the seal.
5. Carry out the procedure in reverse to reassemble the components.
6. Check that the spray operates correctly and that there are no leaks.


Fig. 8 Reassembly of Filter

Warning

Replace a worn or damaged filter cartridge.

7. Troubleshooting

If the device does not operate correctly review the checklist below in order to exclude any user error before contacting your local supplier or VDW GmbH directly.

Problem	Possible cause	Solution
No operation	Defective power cord connection	– Check mains socket – Return to VDW After Sales Service
	Blown fuse	Return to VDW After Sales Service
No water spray	Defective water nose connection	Check water supply
	No water pressure	Check pressure of mains
	Clogged filter	Clean or replace filter
	Defective solenoid valve	Return to VDW After Sales Service
No water at tip but vibration present	Obstructed tip or irrigation file	Remove obstruction for tips or irrigation files
	Wrong tip selection	Check tip
	Incorrect spray adjustment	Re-adjust water spray
Insufficient power	Worn or deformed tip	Replace tip
Low vibration	Presence of fluid or moisture between handpiece and cord	Thoroughly dry electrical contacts with air
No ultrasonic output	Incorrect tip tightening	Re-tighten with wrench
	Defective connector contact	Clean cord contacts
	Cut wire in handpiece cord	Return to VDW After Sales Service for cord replacement
Water leak between handpiece and its base or at junction between handpiece or cord	Worn O-ring in handpiece	Replace the O-ring

Note

Find additional information in the *VDW Endo Tip Guide* giving specific recommendations on how to properly use the *VDW endodontic tips*.

8. Technical / Performance Data

Model/Type	VDW.ULTRA® Ultrasonic device
Minimum frequency of vibrations	28 KHz
Weights	Control unit: 0.9 kg Handpiece: 52 gr
Dimensions (W x D x H)	Control unit: 11.2cm x 11.5cm x 12cm Foot pedal: 7cm x 9.5cm x 3cm
Maximum power	30 VA
Mains power supply	100 VAC / 115 VAC / 220 VAC / 230 VAC 50 Hz / 60 Hz
Protection	Voltage of 100 VAC and 115 VAC: values of fuses 1 and 2 are 500 mA (not accessible) Voltage of 220 VAC and 230 VAC: values of fuses 1 and 2 are 315 mA (not accessible) Values of fuse 3 is 1.5 AT (not accessible)
Device protection class	Protection Class II
Applied component type	BF
Intermittent operation	5 mins. Operation / 10 mins. Pause
Water pressure at inlet	1 to 5 bar (14 to 72 PSI)
Water filter	60 µm
Turn half way to adjust water flow	0-90 ml/min
Ambient conditions	
Operating temperature	+10 to +40°C at a maximum of 95% relative humidity without condensation
Storage temperature	-20 to +70°C at a maximum of 95% relative humidity without condensation

Note

The voltage is preset. The device is protected by fuses (unaccessible for user).

9. Disposal and Recycling

As VDW.ULTRA® falls under the category of electric/electronic equipment, the unit must be disposed of according to a specialized procedure for collection, pick-up and recycling or destruction (especially on the European market, with reference to the Directive 2002/96/EEC of 27.01.2003).

As soon as the unit reaches the end of its life cycle, it is strongly recommended to contact the local dealer to obtain information on how to proceed.

10. Warranty

Important – please read this limited warranty carefully to understand your rights and obligations.

The manufacturer hereby certifies the correct construction of the product, use of top-quality materials, performance of all the required testing, and product compliance to pertinent laws and regulations in force.

The product is covered by warranty for a period of **12 months (excluding tips)** beginning on the date of delivery to the end user as attested by the billing documents, containing the serial number of the unit, issued by the seller at the moment of purchase.

The user has the right to request assistance under warranty only during the warranty period and only on condition that he notifies the manufacturer of the defect in writing within 2 months from the date the defect is noticed.

1. This warranty covers only replacement or repair of single components or pieces affected by manufacturing defects, excluding all expenses for visits by technical personnel and packing expenses.

2. This warranty does not extend to damage or breakdowns attributable to incorrect maintenance, incorrect power supply, negligence, and/or unskilled handling. The warranty does not cover trouble attributable to the user's negligent failure to perform ordinary maintenance (see instruction manual).

3. This warranty does not provide for any compensation for direct or indirect damages of any kind to persons and/or things attributable to the inefficient operation of the unit.

4. This warranty shall automatically be void should the unit be repaired, modified, or in any manner tampered with by the user or by unauthorized third parties.

5. For repairs under warranty, the user must contact only the dealer or the manufacturer's authorized service centres, or the manufacturer themselves. It is the responsibility of the user to ensure that the packing is in perfect condition or in any case suitable for safely shipping the unit.

6. Controversies regarding application of the warranty or the quality or condition of the devices supplied shall not grant the purchaser the faculty to suspend payment for the devices, in full or in part by installment.

7. The purchaser shall have no right to request compensation for device down-time.

8. This warranty expressly excludes defects:

- caused by damage during shipping.

- which are not attributable to manufacturing defects but rather to normal wear caused by proper use of the device (and in particular, calcification and wear of the parts subject to wear).

- caused by atmospheric events such as (for example) lightning, by fires, and by humidity and in general by any other situation or event not directly attributable to the manufacturer.

9. This warranty does not include transportation costs.

The warranty has validity only if accompanied by the billing document proving date of purchase of the product.

Note

The manufacturer has the right to modify the device and/or the user manual without prior notice.

Service Form

In case your device is defective please return it to VDW and include this Service Form.

Your address (or stamp):

Name: _____

Street: _____

Post code, place, country: _____

Telephone: _____

Fax: _____

E-mail: _____

Contact person: _____

The device was purchased on (date): _____

at (e.g. dealer, or sales rep.): _____

Please clean and sterilize the device before returning it, and include all accessories:

VDW GmbH – After Sales Service

Bayerwaldstraße 15, 81737 München

Description of defect: _____

- Repair within the warranty (invoice necessary)
- If repair costs exceed _____ € please send an estimate.
- Please repair the device, I do not need an estimate.

Date: _____ Signature: _____

APPENDIX

11. Electromagnetic Compatibility

Notes

- The VDW.ULTRA® requires special precautions with regard to electromagnetic compatibility.
- It must be installed and prepared for use as described in section 4 "Operation".
- Certain types of mobile telecommunication devices such as mobile telephones are likely to interfere with the VDW.ULTRA®.
- The recommended separation distances in this paragraph must therefore be complied with.
- The VDW.ULTRA® must not be used near or on top of another device. If this cannot be avoided, it is necessary – before clinical use – to check the equipment for correct operation under the conditions of use.
- The use of accessories other than those specified or sold by VDW as replacement parts may have the consequence of increasing the emissions or decreasing the immunity of the unit.

11.1. Electromagnetic Emissions

Notes

- VDW.ULTRA® is intended for use in the electromagnetic environment specified in the following tables 1, 2 and 3 below.
- The user and/or installer of the unit must ensure that it is used in such an environment.

Table 1

Emission test		Electromagnetic environment - guidance
RF emission - CISPR 11	Group 1	VDW.ULTRA® uses RF energy for internal operation. Therefore, its radiofrequency emissions are very low and are not likely to cause any interference in nearby equipment.
RF emission- CISPR 11	Class B	
Harmonic current emission IEC 61000-3-2	Class A	VDW.ULTRA® is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes
Voltage fluctuation and flicker IEC 61000-3-3	Complies	

11.2. Electromagnetic Immunity

Table 2

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 6\text{KV}$ contact $\pm 8\text{KV}$ air	$\pm 6\text{KV}$ contact $\pm 8\text{KV}$ air	Floors must be wood, concrete, cement or tiled. If floors are covered with synthetic material (carpet, etc.), the relative humidity must be at least 30%.
Electrical fast transients IEC 61000-4-4	$\pm 2\text{KV}$ for power supply lines	$\pm 2\text{KV}$ for power supply lines	Power quality should be that of a typical commercial or hospital environment (hospital, clinic).
Surges IEC 61000-4-5	$\pm 1\text{KV}$ differential mode $\pm 2\text{KV}$ common mode	$\pm 1\text{KV}$ differential mode $\pm 2\text{KV}$ common mode	Power quality should be that of a typical commercial or hospital environment (hospital, clinic).
Voltage dips, short interruptions and voltage variations IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 250 cycles	<5% UT (>95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 250 cycles	Power quality should be that of a typical commercial or hospital environment. If the use of the VDW.ULTRA® requires continued operation during a power cut, it is recommended that the product be powered from a separate power supply (UPS, etc.).

11.3. Electromagnetic Immunity / Mobile Radiofrequency Equipment

Table 3

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Portable and mobile radiofrequency communications devices must not be used near the VDW.ULTRA® (including its cables) at a distance less than that recommended and calculated according to the frequency and power of the emitter.</p>			
Conducted disturbance, radiofrequency fields. IEC 61000-4-6	3 V/m 150 KHz to 80 MHz	3 V/m	Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated radiofrequency electromagnetic field. IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum power rating of the emitter in watts (W) according to the manufacturer's specifications and d is the recommended minimum separation distance in meters (m).
<p>The electromagnetic field strengths of fixed radiofrequency emitters, as determined by an electromagnetic environment measurement (a), must be less than the compliance level in each frequency range (b).</p> <p> Interference may occur near equipment marked with this symbol.</p> <p> Note: At 80 MHz and 800 MHz, the higher frequency range applies.</p>			

 **Note**

These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

– The electromagnetic field strengths of fixed radiofrequency emitters, such as base stations for mobile telephones (cellular/cordless), mobile radios, amateur radios, AM/FM radio broadcasts and TV broadcasts cannot be determined exactly by theory. To assess the electromagnetic environment due to fixed radiofrequency emitters, an electromagnetic environment measurement must be made. If the measured radiofrequency field strength in the immediate environment where the product is used exceeds the compliance level specified above, the performance of the product must be tested to verify whether it conforms to the specifications. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product.

– In the 150 kHz to 80 MHz frequency range, the electromagnetic field strengths must be less than 3 V/m.

11.4. Recommended Separation Distances

The VDW.ULTRA® is intended for use in an electromagnetic environment in which radiated radiofrequency disturbances are controlled.

The user and/or installer of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radiofrequency communications equipment (emitters) and the VDW.ULTRA®, according to the maximum output power of the equipment, as recommended in the table below.

Table 4

Rated max. power of the emitter (W)	Separation distance in meters (m) according to emitter frequency		
	150 KHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0.01	0.12 m	0.12 m	0.23 m
0,1	0.38 m	0.38 m	0.73 m
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m

For emitters rated at max. power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the emitter, where P is the max. power rating of the emitter in watts (W) according the manufacturer.

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note

These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

11.5. Cable Lengths

Both handpiece cord and power cord are of 2m length.

Table 5

Cables and accessories	Maximum length	Complies with
Handpiece cord Control foot pedal Power cord	<3m	RF emission, CISPR 1 – Class B/Group 1 Harmonic current emission: IEC 61000-3-2 Voltage fluctuation: IEC 61000-3-3 Immunity to electrostatic discharge: IEC 61000-4-2 Immunity to electrical fast transients/bursts: IEC 61000-4-4 Immunity to surges: IEC 61000-4-5 Immunity to voltage dips, short interruptions and voltage variations: IEC 61000-4-11 Immunity to conducted disturbances induced by radiofrequency fields: IEC 61000-4-6 Immunity to radiated radiofrequency electromagnetic fields: IEC 61000-4-3

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