

## Essential endo motor

### User-friendly

- Torque controlled motor for safe operation
- Broad torque range for operation of several NiTi systems
- Individual selection of torque and speed values for all NiTi systems
- Additional individual program for storage of max. 15 personal settings
- Convenient battery operation, motor can also be used while charging



## → Focus on treatment

- Pre-programmed torque and speed settings for Mtwo®, FlexMaster® and Gates
- Intuitive menu navigation
- Simple operation and clear display



## → Automatic safety functions

- Automatic reverse rotation when set torque limit is reached
- Acoustic warning signals
  - at 75 % of the set torque values
  - when rotating in reverse direction

## → High quality contra-angle

- Small head for good visibility and easy handling
- Fast file change



### Range of products

#### REF 1153 VDW.SILVER®

Including Sirona 6:1 contra-angle and foot control

#### REF 1079 Separate Sirona 6:1 contra-angle with small head

#### REF 1155 VDW.SILVER® (REF 1153) + Mtwo® StarterKit

Contents Mtwo® StarterKit:  
2 Mtwo® of sizes 10/.04, 15/.05, 20/.06, 25/.06 each  
1 Mtwo® of sizes 30/.05, 35/.04, 40/.04, 25/.07 each  
1 Mtwo® SystemBox  
3 plastic training blocks  
2 x 3ml FileCare® EDTA  
1 user CD-ROM and instruction brochure

#### REF 1154 VDW.SILVER® (REF 1153) + FlexMaster® StarterKit

Contents FlexMaster® StarterKit:  
Each 1 FlexMaster® BasicKit 21 and 25 mm  
1 FlexMaster® SystemBox  
3 plastic training blocks  
2 x 3ml FileCare® EDTA  
1 instruction brochure